exaSHAPE

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




Bioshield Pocket® is also available in a **KIT version**, including 2 resorbable monofilament PGCL sutures for the preparation of the "purse string" assembly.

REF	DESCRIPTION	SIZE (w x h x th)
POCKET S 4.0	KIT Pocket small	18 x 14 cm x 0,5 mm
POCKET M 4.0	KIT Pocket medium	20 x 16 cm x 0,5 mm
POCKET L 4.0	KIT Pocket large	22 x 18 cm x 0,5 mm
POCKET XL 4.0	KIT Pocket extra large	23 x 20 cm x 0,5 mm



Bioripar® - Exashape PREPEC

REF	DESCRIPTION	SIZE (w x h x th)
AEPB(F)058-100S	Rectangular, meshed	10 x 5 cm x 0,5 mm
AEPB(F)058-200S	Rectangular, meshed	20 x 5 cm x 0,5 mm
AEPB(F)108-150S	Rectangular, meshed	15 x 10 cm x 0,5 mm
AEPB(F)108-200S	Rectangular, meshed	20 x 10 cm x 0,5 mm
AEPB(F)158-150S	Squared, meshed	15 x 15 cm x 0,5 mm
AEPB(F)158-200S	Rectangular, meshed	20 x 15 cm x 0,5 mm

Prepectoral

Portafolio



Bioripar® - Exashape GRID

REF	DESCRIPTION	SIZE (w x h x th)
AEPB(F)084-162S	Half-moon, perforated	16 x 8 cm x 0,5 mm
AEPB(F)074-172S	Half-moon, perforated	17 x 7 cm x 0,5 mm
AEPB(F)104-202S	Half-moon, perforated	20 x 10 cm x 0,5 mm

REF	DESCRIPTION	SIZE (w x h x th)
AEPB(F)104-150S	Rectangular, perforated	15x 10 cm x 0,5 mm
AEPB(F)104-200S	Rectangular, perforated	20x 10 cm x 0,5 mm

Bioripar® - Exashape EXPANDER


REF	DESCRIPTION	SIZE (w x h x th)
AEPB(F)064-080S	Rectangular, perforated	8 x 6 cm x 0,5 mm

Portafolio


- PROTECTION OF UPPER QUADRANTS
- PROTECTION OF LOWER QUADRANTS
- BIOSHIELD PRO IN CONJUNCTION WITH BIOSHIELD POCKET

• Rehydrate the Exashape membrane in sterile saline solution.

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



Welcome to the 3rd Dimension

ULTRA LIGHT

Less than 0.6 mm thick, this membrane helps to minimize the foreign body response by facilitating rapid integration and faster tissue regeneration. The overall biological mass is 50% less than other membranes dermal based.

PRE-SHAPED

The membrane's shape is conceived for anatomical and round implants. A perfect support is achieved by tightening the petals. wrinkle-free, tailor-made protection for each implant profile or volume, without waving or wrinkling which, as with thick matrices, can be sometimes perceptible in thinner patients

MESHED

The meshed scaffold design optimizes fluid drainage preventing accumulation. The perfect balance in the proportion between cut-outs and collagen best support timing of regeneration

FAST, TOUCHLESS ASSEMBLY

The new pre-shape design facilitates a rapid, touchless assembly procedure that can be completed in less than 1 minute minimizing the risk of contamination and ensuring a sterile environment for optimal surgical outcomes

SELECTIVE SHIELD

Bioshield pocket 3D adapts more effectively to the implant and it shields only where needed: the dermal flap interface. This way foreign body reaction is minimized, and also healing process is not overloaded with unnecessary effort to remodeling exceeding biomass.

REFERENCES:

- 1 Bielli, A., Bernardini, R., et al. (2018) Characterization of a new decellularized bovine pericardial biological mesh: Structural and mechanical properties. *Journal of the Mechanical Behavior of Biomedical Materials*. 78 (2018) 420-426
- 2 Bernardini, R., Varvaras D., D'Amico F., et al. (2019) Biological acellular pericardial mesh regulated tissue integration and remodeling in a rat model of breast prosthetic implantation. *J Biomed Mater Res*. 2019:1-14.
- 3 Varvaras, D., et al. (2017) Safety, tolerability, and efficacy evaluation of immediate total wrapping with biological mesh implant-based breast reconstruction: an under-estimated subcutaneous approach with "biological texturization" prostheses. *Preclinical animal study. The Gulf Journal of Oncology, Supplement January 2017*.
- 4 Capuano, I., Bernardini, R., Varvaras, D., Mattei, M. (2020) Acellular Dermal Matrix in Prosthetic Breast Reconstructive Surgery with Prepectoral Technique: A Literature Review. *Journal of Experimental Pathology*. Volume 1, Issue 2: 50-59

exASHAPE

Bioshield
Pocket
3D

PRE-SHAPED
READY IN 60"
RAPIDLY INTEGRATED

REF

DESCRIPTION

SIZE (W x H x Proj.)

AEPB(F)154-158H050	Pocket 3D Size A	15,5 x 15 x 5 cm
AEPB(F)174-178H055	Pocket 3D Size B	17,5 x 17 x 5,5 cm
AEPB(F)194-198H065	Pocket 3D Size C	19,5 x 19 x 6,5 cm

THIS CATALOGUE IS INTENDED FOR HEALTHCARE PROFESSIONALS.
Refer to the Instructions for Use for detailed information on Intended Use, Warnings and Precautions.

INTERNATIONAL DISTRIBUTOR
Advanced Biomedical Concept Srl
Via Sabotino 2, 00195 Roma, Italia
www.advancedbioconcept.it
+39 06 86357956

MANUFACTURER
 Assut Europe SpA
Via Giuseppe Gregoraci, 12
00173 Roma, Italia
www.assuteurope.com

2025

BioRipar®

exaSHAPE

Estudios

Estudios *in vitro*

Characterization of a new decellularized bovine pericardial biological mesh: Structural and mechanical properties

Alessandra Bielli^{a,1}, Roberta Bernardini^{b,e,1}, Dimitrios Varvaras^{c,1}, Piero Rossi^c,
Giancarlo Di Blasi^d, Giuseppe Petrella^c, Oreste Claudio Buonomo^c, Maurizio Mattei^{a,e,*,1},
Augusto Orlandi^{a,1}

^a Institute of Anatomic Pathology, Department of Biomedicine and Prevention, Tor Vergata University, Rome, Italy

^b Centro Servizi Interdipartimentale – STA, University of Rome “Tor Vergata”, Rome, Italy

^c Department of Experimental Medicine and Surgery, University of Rome “Tor Vergata”, Italy

^d Brachi Testing Services Srl, Prato, Italy

^e Dept. of Biology, University of Rome “Tor Vergata”, Rome, Italy

Table 1

Comparison of thickness and mechanical function of decellularized bovine pericardium.

Mesh	Tissue Thickness (mm)	Tensile Stress (N)	Extension (%)
Bioripar* (n = 4)	0.50 ± 0.04	22.9 ± 6.8*	25.7 ± 6.1
Tutomesh* (n = 5)	0.50 ± 0.05	7.7 ± 3.9	18.5 ± 5.6

n = number of specimens; * = $p < 0.05$.

Table 2

Bursting strength test.

	Bursting strength (KPa)
Bioripar* (n = 3)	560 ± 0*
Tutomesh* (n = 3)	180 ± 19

n = number of specimens for mesh.

* Bioripar* did not burst at the maximal scale value.

- 2018
- **Propiedades mecánicas**
Exashape vs Tuttomesh.
- Exashape muestra tener una **mayor resistencia** que Tuttomesh

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

Biological acellular pericardial mesh regulated tissue integration and remodeling in a rat model of breast prosthetic implantation

Roberta Bernardini^{1†} | Dimitrios Varvaras^{2†} | Federico D'Amico^{3†} |
 Alessandra Bielli^{3†} | Maria Giovanna Sciolli³ | Filadelfo Coniglione^{4,5} | Piero Rossi² |
 Oreste C. Buonomo² | Giuseppe Petrella² | Maurizio Mattei^{1,6†} | Augusto Orlandi^{3,5†}

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FOR INTERNAL USE ONLY

Estudios clínicos

Breast Center
Villa Tiberia Hospital
Roma

Pre-Pectoral Breast Reconstruction And Complete Implant Coverage With A New Bovine Acellular Pericardium Matrix: Breast Center Villa Tiberia Hospital Experience

*D. Varvaras^{*1}, M.G. Lo Schiavo², V. De Felice¹, S. Tardioli³, G. Papadopoulos⁴,
G. Trobbiani⁵, R.A. Madaio⁶, H. Nina⁷, T.M.J. Al Lawati⁸.*

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THE USE OF EXASHAPE™ BIOSHIELD POCKET IN PREPECTORAL BREAST RECONSTRUCTION: A PRELIMINARY EXPERIENCE

Marco Mazzocchi, Giorgio Sigorini, Eugenio Cerciello,
Francesco Barberini

Breast Unit, University of Perugia, Perugia, Italy

Figure 2. A bilateral mastectomy case: pre-operative view (on the left) and post-operative view (on the right) two months after surgery.

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Table II. Surgical characteristic of the patients enrolled in the study.

Unilateral mastectomy	5
Bilateral mastectomy	16
Skin sparing mastectomy	15
Nipple sparing mastectomy	6
Exashape® size	12 medium size and 25 large size
Prosthesis volume (mL)	200-500 (mean value: 315)
Day with drainage	7-15 (mean value: 10)

Skin-reducing mastectomy and immediate hybrid breast reconstruction using combination of acellular pericardium matrix and de-epithelialized dermal flap in patients with medium and large hypertrophic and ptotic breasts

D. Varvaras¹, A. Sachoulidou², V. De Felice¹, F. Frusone¹, G. Papadopoulos³, S. Tardioli⁴, S. D'Onofrio⁵, T. Al Lawati⁶, H. Nina⁷

¹Surgery, Tiberia Hospital, Rome, Italy, ²Second Department of Propaedeutic Surgery, Aristotle University of Thessaloniki, Thessaloniki, Greece, ³Breast Imaging Division - Department of Radiology, Tiberia Hospital, Rome, Italy, ⁴Breast Imaging Division - Department of Radiology, Casa di Cura INI, Grottaferrata, Italy, ⁵Euromedica Geniki Kliniki, Thessaloniki, Greece, ⁶Surgery, Royal Hospital, Muscat, Oman, ⁷Breast Imaging Division - Department of Radiology, San Carlo di Nancy, Rome, Italy, ⁸Surgery, University Medical Center of Tirana Mother Teresa, Tirana, Albania

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

Use of Acellular Pericardial Biological Mesh for Prepectoral and Dual Plane Prosthetic Breast Reconstruction Following Conservative Mastectomies

D. Varvaras¹, F. Frusone¹, A. Sachoulidou², S. Tardioli³, S. D'Onofrio⁴, G. Papadopoulos⁵, T. Al Lawati⁶, V. De Felice¹, L. La Barbera⁷, H. Nina⁸

¹Surgery, Tiberia Hospital, Rome, Italy, ²Second Department of Propaedeutic Surgery, Aristotle University of Thessaloniki, Thessaloniki, Greece, ³Breast Imaging Division - Department of Radiology, Tiberia Hospital, Rome, Italy, ⁴Breast Imaging Division - Department of Radiology, Casa di Cura INI, Grottaferrata, Italy, ⁵Euromedica Geniki Kliniki, Thessaloniki, Greece, ⁶Surgery, Royal Hospital, Muscat, Oman, ⁷Breast Imaging Division - Department of Radiology, San Carlo di Nancy, Rome, Italy, ⁸Surgery, University Medical Center of Tirana Mother Teresa, Tirana, Albania

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A Pericardium Bovine Matrix Pocket in DTI Prepectoral Breast Reconstruction

Roy De Vita,¹ Maria Lucia Mangialardi,^{1,a} Marcello Pozzi,¹ Veronica Vietti
Michelina,¹ Jacopo Nanni,¹ Nicola Zingaretti,² Pier Camillo Parodi,²
Giovanni Zoccali¹

Figure 2 Patient with previous QUART of the right breast, submitted to prophylactic bilateral NSM for BRCA2 mutation + prepectoral implant ABPM assisted immediate reconstruction (anatomical microtextured implants, moderate projection; 415 cc). Preoperative (first row) and 1 year postoperative (second row).


Figure 3 Patient submitted to prophylactic bilateral NSM for BRCA2 mutation + prepectoral implant ABPM assisted immediate reconstruction (anatomical microtextured implants, moderate projection; 420 cc). Preoperative (first row) and 1 year postoperative (second row).

A Pericardium Bovine Matrix Pocket in DTI Prepectoral Breast Reconstruction

Roy De Vita,¹ Maria Lucia Mangialardi,^{1,a} Marcello Pozzi,¹ Veronica Vietti Michelina,¹ Jacopo Nanni,¹ Nicola Zingaretti,² Pier Camillo Parodi,² Giovanni Zoccali¹

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Table 5 Current Available Data on Complications Rate Related to the Use of ABPM in Breast Reconstruction

Study	n°pts (breasts)	Implant Positioning	Seroma	Wound Deishence	Infection	Skin Necrosis	Implant Loss	Capsular Contracture
Borgognone ²⁵ 2011	1 (2)	RP	-	-	-	-	-	-
Semprini ²⁶ 2012	1 (2)	RP	-	-	-	-	-	-
Modir ¹⁴ 2012	54 (93)	RP	7	-	6	5	2	-
Dawson ²⁷ 2013	7 (12)	RP	-	-	-	-	-	-
Gubitos ¹⁵ 2014	24 (28)	RP	5	-	2	-	1	-
Eichler ¹⁶ 2017	45 (27) Comparative with 18 cases ADM	RP	-	-	5	0	1	6
Castagnetti ¹⁷ 2020	123 (141)	RP	4	5	2	5	3	1
Wang ¹⁸ 2021	100 (44) Comparative with 45 cases LD assisted; 11, 2 stages	RP	-	-	1	3	2	1
Mazzocchi ¹⁹ 2022	21 (37)	PP	1	4	0	0	1	0

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

BREAST

Complications of Acellular Dermal Matrices in Breast Surgery

Ron Israeli, M.D.
Hempstead, N.Y.

Summary: Acellular dermal matrices have been used in breast surgery for a decade. They are widely used in implant-based breast reconstruction to provide coverage of the inferolateral aspects of the prosthesis. Numerous benefits have been reported with this approach including improved fold control, better support and control of the implant pocket with concomitant reduced risk of malposition, and improved lower pole expansion. Seroma, infection, mastectomy skin necrosis, and expander/implant loss are the most commonly reported complications with this approach, and the incidences vary widely among studies. Patient selection and adherence to established intraoperative technique principles related to acellular dermal matrix use are both critical to minimizing the risk of complications. Acellular dermal matrices are also being used in aesthetic breast surgery, revision breast surgery, and nipple reconstruction, but clinical experience is limited. This article reviews the complications associated with the use of matrices in breast surgery from the published literature. (*Plast. Reconstr. Surg.* 130 (Suppl. 2): 159S, 2012.)

The first reported use of acellular dermal matrix in breast surgery was in 2001 for the correction of implant rippling.¹ Over the past decade, the indication for the use of acellular dermal matrices in breast surgery has expanded to include primary implant-based reconstruction, aesthetic breast surgery, revision breast surgery, and nipple reconstruction. Their use is now common practice in implant-based breast reconstruction.² Despite their popularity, an ongoing concern is the risk of complications, particularly seroma, skin necrosis, infection, and failed reconstruction with expander or implant loss. An increased rate of these complications has been reported in some recent studies that have questioned the net benefit of acellular dermal matrix use.³⁻⁵

This article reviews the incidences and types of complications associated with the use of acellular dermal matrices for the various indications in breast surgery. Risk factors that may influence complications as well as strategies that may help reduce the risk of complications are also discussed.

From the Department of Surgery, Division of Plastic Surgery, Hofstra University School of Medicine in Partnership with the North Shore-LIJ Health System.
Received for publication February 1, 2012; accepted March 26, 2012.
Copyright © 2012 by the American Society of Plastic Surgeons
DOI: 10.1097/PRS.0b013e3182634e62

www.PRSJournal.com

159S

BREAST

Embrace the Change: Incorporating Single-Stage Implant Breast Reconstruction into Your Practice

Jose Rodriguez-Felix, M.D.
Mark A. Codner, M.D.

Atlanta, Ga.

Background: Multiple studies have reported on the safety of nipple-sparing mastectomy and low complication rates associated with single-stage implant breast reconstruction. Yet many plastic surgeons continue to be resistant to change. This article presents the senior author's (M.A.C.) experience during his transition period from the latissimus dorsi flap with adjustable implants to a "one-and-done" approach using shaped implants and fetal bovine acellular dermal matrix.

Methods: A literature review was performed selecting articles discussing single-stage implant reconstruction, indications, outcomes, technique, and complications. Additional articles were selected after review of the references of identified articles. Clinical pearls discussed include patient selection, implant selection, and mastectomy incision choices, with a detailed description of the senior author's operative technique.

Results: Twenty-seven single-stage implant reconstructions were performed. Average mastectomy weight was 343.82 g. The average implant volume was 367 cc. Shaped implants were most commonly used. Acellular dermal matrix was used in all breasts. Complications included erythema requiring intravenous antibiotics (three patients), skin ischemia caused by methylene blue (one patient), seroma (one patient), unilateral partial nipple necrosis (one patient), mastectomy skin necrosis (one patient), and exposed/infectious implants that were salvaged using a sequential irrigation protocol described by Sforza et al. in 2014 (two patients).

Conclusions: Breast reconstruction after mastectomy has evolved toward less invasive, single-stage procedures. Aesthetic refinements include nipple-sparing mastectomy, use of acellular dermal matrix, shaped implants, and fat grafting. Selected patients will benefit from a one-and-done breast implant reconstruction with no additional oncologic risk. Surgeons must embrace the change and provide their patients with a procedure that will offer the best aesthetic outcomes. (*Plast. Reconstr. Surg.* 136: 221, 2015.)

CLINICAL QUESTION/LEVEL OF EVIDENCE: Therapeutic, IV.

One in eight women in the United States will develop invasive breast cancer during their lifetime. An estimated total of 292,670 women will be diagnosed this year (2014) with invasive breast cancer and another 62,570 with carcinoma in situ.¹ Approximately 40 percent of women diagnosed with breast cancer will be treated with mastectomy, with fewer than 25 percent undergoing immediate reconstruction.²

From Emory University and Mark Codner MD Plastic Surgery.
Received for publication December 14, 2014; accepted February 10, 2015.
Copyright © 2015 by the American Society of Plastic Surgeons
DOI: 10.1097/PRS.0000000000001448

Disclosure: Dr. Rodriguez-Felix has no financial interest in any of the products or devices mentioned in this article. Dr. Codner is a consultant speaker for Mentor Corporation and Ulthera, Inc., and receives royalties from publishers QMP and Elsevier. He receives a fellowship grant from TEI Biosciences, Inc.

Supplemental digital content is available for this article. Direct URL citations appear in the text; simply type the URL address into any Web browser to access this content. Clickable links to the material are provided in the HTML text of this article on the Journal's website (www.PRSJournal.com).

www.PRSJournal.com

221

Hindawi Publishing Corporation
Plastic Surgery International
Volume 2014, Article ID 172004, 10 pages
http://dx.doi.org/10.1155/2014/172004

Review Article

Acellular Dermal Matrices and Radiotherapy in Breast Reconstruction: A Systematic Review and Meta-Analysis of the Literature

Luigi Valdatta,^{1,2} Anna Giulia Cattaneo,¹ Igor Pellegatta,² Stefano Scamoni,² Anna Minuti,¹ and Mario Cherubino^{1,2}

¹ Department of Biotechnology & Life Science (DBSV), University of Insubria, 21100 Varese, Italy

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Correspondence should be addressed to Igor Pellegatta; igor.pellegatta@me.com

Received 4 February 2014; Accepted 22 April 2014; Published 21 May 2014

Academic Editor: Georg M. Haemmer

Copyright © 2014 Luigi Valdatta et al. This is an open access article distributed under the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

The increasing use of commercially available acellular dermis matrices for postmastectomy breast reconstruction seems to have simplified the surgical procedure and enhanced the outcome. These materials, generally considered to be highly safe or with only minor contraindications due to the necessary manipulation in preparatory phases, allow an easier one phase surgical procedure, in comparison with autologous flaps, offering a high patient satisfaction. Unfortunately, the claim for a higher rate of complications associated with irradiation at the implant site, especially when the radiation therapy was given before the reconstructive surgery, suggested a careful behaviour when this technique is preferred. However, this hypothesis was never submitted to a crucial test, and data supporting it are often discordant or incomplete. To provide a comprehensive analysis of the field, we searched and systematically reviewed papers published after year 2005 and registered clinical trials. On the basis of a meta-analysis of data, we conclude that the negative effect of the radiotherapy on the breast reconstruction seems to be evident even in the case of acellular dermis matrices aided surgery. However, more trials are needed to make solid conclusions and clarify the poor comprehension of all the factors negatively influencing outcome.

1. Introduction

The acellular dermis matrices (ADM) are products derived from the skin, deprived of their cellular component by standardized treatments [1]. They provide a lower-lateral coverage and support of the implants in the immediate expander/implant-based breast reconstruction after mastectomy. Additional main indications for their use are lack of muscular coverage and cancer invasion to the pectoralis major muscle, and skin nipple sparing mastectomy is a relative indication [2].





The customized commercial products mainly differ in their origin and in procedures for processing, storing, and preparing them before usage. A recent paper compared seven customized ADM suitable for the reconstruction of





the breast, in order to evaluate their cost/benefit ratios, contraindications, and possible side effects [3]. These authors discuss main contraindications for the use of commercial products, among others, namely, the presence of residues of antibiotics or allergenic substances, the lack of sterility, lower strength and elasticity, and higher cost. The Alloderm (LifeCell Corporation, Branchburg, NJ, USA), a customized derivative of the banked human skin, is the most widely used material, despite a few disadvantages: longer rehydration time, possible presence of antibiotics, and nonsterility of the final product.

We analyzed the data reported in peer-reviewed papers, in which irradiation at the site of ADM implant and its timing are considered as possible interfering factors for surgical outcome. The review includes an exhaustive description

Tailoring Braxon[®]


Braxon[®] is pre-shaped in such a way as to contain a breast implant of any size and shape. Its use is intuitive and requires scissors and suture to “dress” the prosthesis and be sutured over the pectoralis major muscle.

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
Table 1
Comparison of thickness and mechanical function of decellularized bovine pericardium.

Mesh	Tissue Thickness (mm)	Tensile Stress (N)	Extension (%)
Bioripar [*] (n = 4)	0.50 ± 0.04	22.9 ± 6.8 [*]	25.7 ± 6.1
Tutomesh [*] (n = 5)	0.50 ± 0.05	7.7 ± 3.9	18.5 ± 5.6

n = number of specimens; ^{*} = p < 0.05.




Table 2
Bursting strength test.

	Bursting strength (KPa)
Bioripar [*] (n = 3)	560 ± 0 [*]
Tutomesh [*] (n = 3)	180 ± 19

n = number of specimens for mesh.

^{*} Bioripar^{*} did not burst at the maximal scale value.

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






Technical Data

- ▶ Material: Titanised type 1a polypropylene mesh
- ▶ Pore size: 1 mm
- ▶ Product weight: 16 g/m²

TiLOOP® Bra Pocket

Recommended fixation points
View of the back of the mesh (which is touching the thoracic wall)

Pre-pectoral position of the TiLOOP® Bra Pocket

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